

Monofocal Plus Intraocular Lens

CE

NMPA



Items	Description
IOL model	PA60AS1
Optic Design	Monofocal, Aspheric
Material	Hydrophobic Acrylic
Lens Design	Single Piece Posterior Chamber Foldable IOL
Haptic Design	L- loop
Haptic Angle	0°
Range of Optical Powers	-10.0D to +40.0D (0.5 D increments)
Refractive Index	1.55
Optic Diameter	6.0mm
Overall Length	13.0mm
A-Constant	(Nominal) A=118.8
	a0=-0.769
	(Haigis) a1=0.234
	a2=0.217
	(HofferQ) pACD=5.64
	(Holl.1) sf=1.84
	(Holl.2) 4.75
	(SRK/T) A=119.0
	(SRK II) A=119.2



Enhanced Intermediate Vision



Improved Depth of Focus



Reduced Glistening

Ventus

Ventus Sağlık Medikal San. Tic. LTD. ŞTİ.

Gazibey Mah. Kızılay Sk. No:4/A Malkara / Tekirdağ

info@ventusmedikal.com www.ventusmedikal.com